





HIP




STAINLESS STEEL HEAD



CoCr HEAD



BIOLOX® FORTE HEAD



BIOLOX® DELTA HEAD

TECHNICAL DOCUMENTATION

12/14
HEADS RANGE




FH ORTHOPEDICS heads range is available in different materials, in several diameters but with a standard 12/14 taper. All the ball heads are compatible with all the FH ORTHOPEDICS stems. Regarding the liners, the ceramic liners can be match only with the ceramic heads either BIOLOX® FORTE or BIOLOX® DELTA. The TRIANON® highly crosslinked polyethylene liner can be match only with CoCr heads or with the BIOLOX® ceramic heads.

REFERENCES - NOT AVAILABLE FOR USA




STAINLESS STEEL HEADS	
231402	STAINLESS STEEL HEAD TAPER 12/14 D.22,2 M (0)
231403	STAINLESS STEEL HEAD TAPER 12/14 D.22,2 L (+2)
247852	STAINLESS STEEL HEAD TAPER 12/14 D.28 XS (-5)
231391	STAINLESS STEEL HEAD TAPER 12/14 D.28 S (-3,5)
231392	STAINLESS STEEL HEAD TAPER 12/14 D.28 M (0)
231393	STAINLESS STEEL HEAD TAPER 12/14 D.28 L (+3,5)
231394	STAINLESS STEEL HEAD TAPER 12/14 D.28 XL (+7)
231386	STAINLESS STEEL HEAD TAPER 12/14 D.32 S (-3,5)
231387	STAINLESS STEEL HEAD TAPER 12/14 D.32 M (0)
231388	STAINLESS STEEL HEAD TAPER 12/14 D.32 L (+3,5)



BIOLOX® FORTE HEADS	
234131	BIOLOX® FORTE HEAD TAPER 12/14 D.28 S (-3,5)
234132	BIOLOX® FORTE HEAD TAPER 12/14 D.28 M (0)
234133	BIOLOX® FORTE HEAD TAPER 12/14 D.28 L (+3,5)
250729	BIOLOX® FORTE HEAD TAPER 12/14 D.32 S (-3,5)
250730	BIOLOX® FORTE HEAD TAPER 12/14 D.32 M (0)
250731	BIOLOX® FORTE HEAD TAPER 12/14 D.32 L (+3,5)
254525	BIOLOX® FORTE HEAD TAPER 12/14 D.36 S (-3,5)
254526	BIOLOX® FORTE HEAD TAPER 12/14 D.36 M (0)
254527	BIOLOX® FORTE HEAD TAPER 12/14 D.36 L (+3,5)



COBALT CHROMIUM HEADS	
240858	COCR HEAD TAPER 12/14 D.22,2 M (0)
240859	COCR HEAD TAPER 12/14 D.22,2 L (+2)
240862	COCR HEAD TAPER 12/14 D.28 S (-3,5)
240863	COCR HEAD TAPER 12/14 D.28 M (0)
240864	COCR HEAD TAPER 12/14 D.28 L (+3,5)
240865	COCR HEAD TAPER 12/14 D.28 XL (+7)
240870	COCR HEAD TAPER 12/14 D.32 S (-3,5)
240871	COCR HEAD TAPER 12/14 D.32 M (0)
240872	COCR HEAD TAPER 12/14 D.32 L (+3,5)
240873	COCR HEAD TAPER 12/14 D.32 XL (+7)
263684	COCR HEAD TAPER 12/14 D.36 S (-3,5)
263685	COCR HEAD TAPER 12/14 D.36 M (0)
263686	COCR HEAD TAPER 12/14 D.36 L (+3,5)
263687	COCR HEAD TAPER 12/14 D.36 XL (+7)



BIOLOX® DELTA HEADS	
256516	BIOLOX® DELTA HEAD TAPER 12/14 D.28 S (-3,5)
256517	BIOLOX® DELTA HEAD TAPER 12/14 D.28 M (0)
256518	BIOLOX® DELTA HEAD TAPER 12/14 D.28 L (+3,5)
256519	BIOLOX® DELTA HEAD TAPER 12/14 D.32 S (-3,5)
256520	BIOLOX® DELTA HEAD TAPER 12/14 D.32 M (0)
256521	BIOLOX® DELTA HEAD TAPER 12/14 D.32 L (+3,5)
266371	BIOLOX® DELTA HEAD TAPER 12/14 D.32 XL (+7)
256522	BIOLOX® DELTA HEAD TAPER 12/14 D.36 S (-3,5)
256523	BIOLOX® DELTA HEAD TAPER 12/14 D.36 M (0)
256524	BIOLOX® DELTA HEAD TAPER 12/14 D.36 L (+3,5)
266372	BIOLOX® DELTA HEAD TAPER 12/14 D.36 XL (+7)